




Certificate of Analysis

CP-Buffer-4

AxyGen's ready to use buffers, precipitants, salts, additives, and detergents provide the most complete and economical solution to crystallization of proteins, peptides, small molecules and nucleic acids. They can be used as a rapid Initial screening method or to accelerate the creation of optimization/follow up screens. Reagents are supplied in sterile bottles (100ml/bottle).

Product Name:	1.0 M BisTriPro, pH 8.5
Description:	BisTriPro titrated to pH 5.5 with HCl
Catalog #:	CP-Buffer-4
Chemical Name:	1,3-Bis[tris(hydroxymethyl)methylamino]propane
Molecular Formula:	$\text{CH}_2[\text{CH}_2\text{NHC}(\text{CH}_2\text{OH})_3]_2$
Molecular Weight:	282.33
CAS Number:	64431-96-5
Beilstein Registry:	1786109
EC Number:	2648993
Purity:	$\geq 99.0\%$
Appearance:	White crystalline powder
Solubility:	Dissolves in water (428 mg/ml), yielding a clear, colorless solution
pKa (25 °C):	(1) 6.8, (2) 9.0
MP:	164-165 °C



AxyGen BioSciences
33210 Central Ave., Union City, CA 94587
Ph.: (510) 476-1090 Fax: (510) 476-1091
Email: sales@axygenbio.com
Website: www.axygenbio.com